




ELECTRO SURGICAL INSTRUMENTS



ALL INSTRUMENTS AVAILABLE IN



MANUFACTURER & EXPORTER



McPHERSON
TC20-035-102
 10.2 cm, 4"

McPHERSON
TC20-038-102
 10.2 cm, 4"
 (Angled)

JEWELER # 5
TC20-041-115
 11.5 cm, 4½"


JEWELER # 5
TC20-044-115
 11.5 cm, 4½"
 (Angled)

JEWELER # 7
TC20-047-115
 11.5 cm, 4½"

ADSON
TC20-050-121
 12.1 cm, 4¾"

| Cat. # | Description | Tip (mm) | Size |
|--------------|----------------------|----------|--------------|
| TC20-035-102 | McPherson | 0.5 | 10.2 cm, 4" |
| TC20-038-102 | McPherson (Angled) | 0.5 | 10.2 cm, 4" |
| TC20-041-115 | Jeweler # 5 | 0.5, 0.7 | 11.5 cm, 4½" |
| TC20-044-115 | Jeweler # 5 (Angled) | 0.5, 0.7 | 11.5 cm, 4½" |
| TC20-047-115 | Jeweler # 7 | 0.5, 0.7 | 11.5 cm, 4½" |
| TC20-050-121 | Adson | 1.0 | 12.1 cm, 4¾" |
| TC20-053-150 | Adson | 1.0 | 15.0 cm, 6" |


For your convenience, you can add the following extensions for tip size at the end of Cat #, Needle Tip (-00), 0.5 mm (-05), 1.0 mm (-10) and 2.0 mm (-20)



BAYONET TC20-056-165 16.5 cm, 6½"
BAYONET TC20-059-180 18.0 cm, 7"
BAYONET TC20-062-200 20.0 cm, 7¾"
BAYONET TC20-071-165 16.5 cm, 6½"
BAYONET TC20-074-180 18.0 cm, 7"
BAYONET TC20-077-200 20.0 cm, 7¾"
BAYONET TC20-086-165 16.5 cm, 6½"
BAYONET TC20-089-180 18.0 cm, 7"
BAYONET 20-092-200 20.0 cm, 7¾"
 (Angled Up) (Angled Up) (Angled Up) (Angled Down) (Angled Down) (Angled Down)

| Cat. # | Description | Tip (mm) | Size |
|--------------|-----------------------|---------------|--------------|
| TC20-056-165 | Bayonet | 0.5, 1.0, 2.0 | 16.5 cm, 6½" |
| TC20-059-180 | Bayonet | 0.5, 1.0, 2.0 | 18.0 cm, 7" |
| TC20-062-200 | Bayonet | 0.5, 1.0, 2.0 | 20.0 cm, 7¾" |
| TC20-065-220 | Bayonet | 0.5, 1.0, 2.0 | 22.0 cm, 8¾" |
| TC20-068-240 | Bayonet | 0.5, 1.0, 2.0 | 24.0 cm, 9½" |
| TC20-071-165 | Bayonet (Angled Up) | 0.5, 1.0, 2.0 | 16.5 cm, 6½" |
| TC20-074-180 | Bayonet (Angled Up) | 0.5, 1.0, 2.0 | 18.0 cm, 7" |
| TC20-077-200 | Bayonet (Angled Up) | 0.5, 1.0, 2.0 | 20.0 cm, 7¾" |
| TC20-080-220 | Bayonet (Angled Up) | 0.5, 1.0, 2.0 | 22.0 cm, 8¾" |
| TC20-083-240 | Bayonet (Angled Up) | 0.5, 1.0, 2.0 | 24.0 cm, 9½" |
| TC20-086-165 | Bayonet (Angled Down) | 0.5, 1.0, 2.0 | 16.5 cm, 6½" |
| TC20-089-180 | Bayonet (Angled Down) | 0.5, 1.0, 2.0 | 18.0 cm, 7" |
| TC20-092-200 | Bayonet (Angled Down) | 0.5, 1.0, 2.0 | 20.0 cm, 7¾" |
| TC20-095-220 | Bayonet (Angled Down) | 0.5, 1.0, 2.0 | 22.0 cm, 8¾" |
| TC20-098-240 | Bayonet (Angled Down) | 0.5, 1.0, 2.0 | 24.0 cm, 9½" |

For your convenience, you can add the following extensions for tip size at the end of Cat #, Needle Tip (-00), 0.5 mm (-05), 1.0 mm (-10) & 2.0 mm (-20)



| | | | | | | | |
|--|---|--|--|--|--|--|---|
| ADSON TC20-101-121 12.1 cm, 4 1/4" | SEMKIN TC20-107-140 14.0 cm, 5 1/2" | CUSHING TC20-113-159 15.9 cm, 6 1/4" | CUSHING TC20-119-178 17.8 cm, 7" | CUSHING (Gursh) TC20-125-178 17.8 cm, 7" | GERALD (Str.) TC20-131-178 17.8 cm, 7" | GERALD (Crv.) TC20-137-178 17.8 cm, 7" | GERALD BAYONET TC20-143-184 18.4 cm, 7 1/4" |
| TC20-104-121 12.1 cm, 4 1/4" | TC20-110-140 14.0 cm, 5 1/2" | TC20-116-159 15.9 cm, 6 1/4" | TC20-122-178 17.8 cm, 7" | TC20-128-178 17.8 cm, 7" | TC20-134-178 17.8 cm, 7" | TC20-140-178 17.8 cm, 7" | TC20-146-184 18.4 cm, 7 1/4" |

| Cat. # | Description | Style | Tip (mm) | Size |
|--------------|-------------|---------------|----------|-----------------|
| TC20-101-121 | Adson | Insulated | 1.0 | 12.1 cm, 4 1/4" |
| TC20-104-121 | Adson | Non-Insulated | 1.0 | 12.1 cm, 4 1/4" |
| TC20-107-140 | Semkin | Insulated | 0.5 | 14.0 cm, 5 1/2" |
| TC20-110-140 | Semkin | Non-Insulated | 0.5 | 14.0 cm, 5 1/2" |
| TC20-113-159 | Cushing | Insulated | 1.0, 1.5 | 15.9 cm, 6 1/4" |
| TC20-116-159 | Cushing | Non-Insulated | 1.0, 1.5 | 15.9 cm, 6 1/4" |
| TC20-119-178 | Cushing | Insulated | 1.0, 2.0 | 17.8 cm, 7" |
| TC20-122-178 | Cushing | Non-Insulated | 1.0, 2.0 | 17.8 cm, 7" |

| Cat. # | Description | Style | Tip (mm) | Size |
|--------------|-----------------|---------------|----------|-----------------|
| TC20-125-178 | Cushing (Gursh) | Insulated | 2.0 | 17.8 cm, 7" |
| TC20-128-178 | Cushing (Gursh) | Non-Insulated | 2.0 | 17.8 cm, 7" |
| TC20-131-178 | Gerald (Str.) | Insulated | 1.0, 1.5 | 17.8 cm, 7" |
| TC20-134-178 | Gerald (Str.) | Non-Insulated | 1.0, 1.5 | 17.8 cm, 7" |
| TC20-137-178 | Gerald (Crv.) | Insulated | 1.0, 1.5 | 17.8 cm, 7" |
| TC20-140-178 | Gerald (Crv.) | Non-Insulated | 1.0, 1.5 | 17.8 cm, 7" |
| TC20-143-184 | Gerald Bayonet | Insulated | 1.0, 1.5 | 18.4 cm, 7 1/4" |
| TC20-146-184 | Gerald Bayonet | Non-Insulated | 1.0, 1.5 | 18.4 cm, 7 1/4" |

For your convenience, you can add the following extensions for tip size at the end of Cat #.
 Needle Tip (-00), 0.5 mm (-05), 1.0 mm (-10) & 2.0 mm (-20)



McPHERSON (Str.)
TC20-149-090
 9.0 cm, 3 1/2"
TC20-152-090
 9.0 cm, 3 1/2"

McPHERSON (Crv.)
TC20-155-090
 9.0 cm, 3 1/2"
TC20-158-090
 9.0 cm, 3 1/2"

COAPTATION NADLER
TC20-161-102
 10.2 cm, 4"
TC20-164-102
 10.2 cm, 4"

COAPTATION
TC20-167-108
 10.8 cm, 4 1/4"
TC20-170-108
 10.8 cm, 4 1/4"

JEWELER # 5
TC20-173-115
 11.5 cm, 4 1/2"
TC20-176-115
 11.5 cm, 4 1/2"

JEWELER # 4
TC20-179-115
 11.5 cm, 4 1/2"
TC20-182-115
 11.5 cm, 4 1/2"

TENZEL
TC20-185-102
 10.2 cm, 4"
TC20-188-102
 10.2 cm, 4"

JEWELER # 7
TC20-191-115
 11.5 cm, 4 1/2"
TC20-194-115
 11.5 cm, 4 1/2"


IRIS (Str.)
TC20-197-108
 10.8 cm, 4 1/4"
20-200-108
 10.8 cm, 4 1/4"



IRIS (Crv.)
20-203-108
 10.8 cm, 4 1/4"
20-206-108
 10.8 cm, 4 1/4"



| Cat. # | Description | Style | Tip (mm) | Size | Cat. # | Description | Style | Tip (mm) | Size |
|--------------|-------------------|---------------|------------|-----------------|--------------|-------------|---------------|------------|-----------------|
| TC20-149-090 | McPherson (Str.) | Insulated | 0.5 | 09.0 cm, 3 1/2" | TC20-179-115 | Jeweler # 4 | Insulated | Micro, 0.5 | 11.5 cm, 4 1/2" |
| TC20-152-090 | McPherson (Str.) | Non-Insulated | 0.5 | 09.0 cm, 3 1/2" | TC20-182-115 | Jeweler # 4 | Non-Insulated | Micro, 0.5 | 11.5 cm, 4 1/2" |
| TC20-155-090 | McPherson (Crv.) | Insulated | 0.5 | 09.0 cm, 3 1/2" | TC20-185-102 | Tenzel | Insulated | Micro, 0.5 | 10.2 cm, 4" |
| TC20-158-090 | McPherson (Crv.) | Non-Insulated | 0.5 | 09.0 cm, 3 1/2" | TC20-188-102 | Tenzel | Non-Insulated | Micro, 0.5 | 10.2 cm, 4" |
| TC20-161-102 | Coaptation Madler | Insulated | 0.5 | 10.2 cm, 4" | TC20-191-115 | Jeweler # 7 | Insulated | Micro, 0.5 | 11.5 cm, 4 1/2" |
| TC20-164-102 | Coaptation Madler | Non-Insulated | 0.5 | 10.2 cm, 4" | TC20-194-115 | Jeweler # 7 | Non-Insulated | Micro, 0.5 | 11.5 cm, 4 1/2" |
| TC20-167-108 | Coaptation | Insulated | 0.5 | 10.8 cm, 4 1/4" | TC20-197-108 | Iris (Str.) | Insulated | 1.0 | 10.8 cm, 4 1/4" |
| TC20-170-108 | Coaptation | Non-Insulated | 0.5 | 10.8 cm, 4 1/4" | TC20-200-108 | Iris (Str.) | Non-Insulated | 1.0 | 10.8 cm, 4 1/4" |
| TC20-173-115 | Jeweler # 5 | Insulated | Micro, 0.5 | 11.5 cm, 4 1/2" | TC20-203-108 | Iris (Crv.) | Insulated | 1.0 | 10.8 cm, 4 1/4" |
| TC20-176-115 | Jeweler # 5 | Non-Insulated | Micro, 0.5 | 11.5 cm, 4 1/2" | TC20-206-108 | Iris (Crv.) | Non-Insulated | 1.0 | 10.8 cm, 4 1/4" |



JANSEN BAYONET
TC20-209-152 15.2 cm, 6"
TC20-212-152 15.2 cm, 6"



JANSEN BAYONET
TC20-215-178 17.8 cm, 7"
TC20-218-178 17.8 cm, 7"
TC20-221-204 20.4 cm, 8"
TC20-224-204 20.4 cm, 8"



MALIS
TC20-227-178 17.8 cm, 7"
TC20-230-178 17.8 cm, 7"



SCOVILLE GREENWOOD
TC20-233-191 19.1 cm, 7½"
TC20-236-191 19.1 cm, 7½"



YASARGIL
TC20-239-197 19.7 cm, 7¾"
TC20-242-197 19.7 cm, 7¾"



HARDY BAYONET
TC20-245-216 21.6 cm, 8½"
TC20-248-216 21.6 cm, 8½"



HARDY BAYONET
TC20-251-216 21.6 cm, 8½"
TC20-254-216 21.6 cm, 8½"
 (Angled)

| Cat. # | Description | Style | Tip (mm) | Size |
|--------------|----------------------|---------------|---------------|--------------|
| TC20-209-152 | Jansen Bayonet | Insulated | 0.5, 1.0, 1.5 | 15.2 cm, 6" |
| TC20-212-152 | Jansen Bayonet | Non-Insulated | 0.5, 1.0, 1.5 | 15.2 cm, 6" |
| TC20-215-178 | Jansen Bayonet | Insulated | 0.5, 1.0, 1.5 | 17.8 cm, 7" |
| TC20-218-178 | Jansen Bayonet | Non-Insulated | 0.5, 1.0, 1.5 | 17.8 cm, 7" |
| TC20-221-204 | Jansen Bayonet | Insulated | 0.5, 1.0, 2.0 | 20.4 cm, 8" |
| TC20-224-204 | Jansen Bayonet | Non-Insulated | 0.5, 1.0, 2.0 | 20.4 cm, 8" |
| TC20-227-178 | Malis | Insulated | 0.5, 1.0, 2.0 | 17.8 cm, 7" |
| TC20-230-178 | Malis | Non-Insulated | 0.5, 1.0, 2.0 | 17.8 cm, 7" |
| TC20-233-191 | Scoville Greenwood | Insulated | 1.0, 1.5 | 19.1 cm, 7½" |
| TC20-236-191 | Scoville Greenwood | Non-Insulated | 1.0, 1.5 | 19.1 cm, 7½" |
| TC20-239-197 | Yasargil | Insulated | 1.0, 1.5 | 19.7 cm, 7¾" |
| TC20-242-197 | Yasargil | Non-Insulated | 1.0, 1.5 | 19.7 cm, 7¾" |
| TC20-245-216 | Hardy Bayonet | Insulated | 1.0 | 21.6 cm, 8½" |
| TC20-248-216 | Hardy Bayonet | Non-Insulated | 1.0 | 21.6 cm, 8½" |
| TC20-251-216 | Hardy Bayonet (Ang.) | Insulated | 1.0 | 21.6 cm, 8½" |
| TC20-254-216 | Hardy Bayonet (Ang.) | Non-Insulated | 1.0 | 21.6 cm, 8½" |

Any European bipolar forceps can be converted into "Irrigation Bipolar Forceps" as shown in previous pages. Add extension "IR" at the end.



BAYONET
TC20-257-01R




BAYONET
TC20-260-01R



BAYONET
TC20-263-01R



BAYONET
TC20-266-01R



ADSON TC20-269-120 12.0 cm, 4 1/4"
McINDOE TC20-272-150 15.0 cm, 6"
GILLIES TC20-275-150 15.0 cm, 6"
DISSECTING TC20-287-150 15.0 cm, 6"
GERALD (1 x 2) TC20-299-178 17.8 cm, 7"
ADSON BAYONET TC20-302-120 16.0 cm, 6 1/4"
ADSON BAYONET TC20-308-160 16.0 cm, 6 1/4" (Ang. Up)
REMINGTON Hobb TC20-317-241 24.1 cm, 9 1/2"

European Monopolar (Male 4.8 mm Connector)

| Cat. # | Description | Tip (mm) | Size |
|--------------|-------------------------|-------------|-----------------|
| TC20-269-120 | Adson | 1.0 mm | 12.0 cm, 4 1/4" |
| TC20-272-150 | McIndoe | 1.0, 2.0 mm | 15.0 cm, 6" |
| TC20-275-150 | Gillies (1 x 2) | -- | 15.0 cm, 6" |
| TC20-278-178 | Gillies (1 x 2) | -- | 17.8 cm, 7" |
| TC20-281-200 | Gillies (1 x 2) | -- | 20.0 cm, 8" |
| TC20-284-240 | Gillies (1 x 2) | -- | 24.0 cm, 9" |
| TC20-287-150 | Dissecting Forceps | 1.0, 2.0 mm | 15.0 cm, 6" |
| TC20-290-178 | Dissecting Forceps | 1.0, 2.0 mm | 17.8 cm, 7" |
| TC20-293-210 | Dissecting Forceps | 1.0, 2.0 mm | 21.0 cm, 8 1/4" |
| TC20-296-241 | Dissecting Forceps | 1.0, 2.0 mm | 24.1 cm, 9 1/2" |
| TC20-299-178 | Gerald (1 x 2) | -- | 17.8 cm, 7" |
| TC20-302-160 | Adson Bayonet | 1.0, 2.0 mm | 16.0 cm, 6 1/4" |
| TC20-305-191 | Adson Bayonet | 1.0, 2.0 mm | 19.1 cm, 7 1/2" |
| TC20-308-160 | Adson Bayonet (Ang. Up) | 1.0, 2.0 mm | 16.0 cm, 6 1/4" |
| TC20-311-191 | Adson Bayonet (Ang. Up) | 1.0, 2.0 mm | 19.1 cm, 7 1/2" |
| TC20-314-210 | Adson Bayonet (Ang. Up) | 1.0, 2.0 mm | 21.0 cm, 8 1/4" |
| TC20-317-241 | Remington Hobb | 1.0, 2.0 mm | 24.1 cm, 9 1/2" |

Monopolar (without Connector)

| Cat. # | Description | Tip (mm) | Size |
|--------------|-------------------------|-------------|-----------------|
| TC20-420-120 | Adson | 1.0 mm | 12.0 cm, 4 1/4" |
| TC20-423-150 | McIndoe | 1.0, 2.0 mm | 15.0 cm, 6" |
| TC20-426-150 | Gillies (1 x 2) | -- | 15.0 cm, 6" |
| TC20-429-178 | Gillies (1 x 2) | -- | 17.8 cm, 7" |
| TC20-432-200 | Gillies (1 x 2) | -- | 20.0 cm, 8" |
| TC20-435-240 | Gillies (1 x 2) | -- | 24.0 cm, 9" |
| TC20-438-150 | Dissecting Forceps | 1.0, 2.0 mm | 15.0 cm, 6" |
| TC20-441-178 | Dissecting Forceps | 1.0, 2.0 mm | 17.8 cm, 7" |
| TC20-444-210 | Dissecting Forceps | 1.0, 2.0 mm | 21.0 cm, 8 1/4" |
| TC20-447-241 | Dissecting Forceps | 1.0, 2.0 mm | 24.1 cm, 9 1/2" |
| TC20-450-178 | Gerald (1 x 2) | -- | 17.8 cm, 7" |
| TC20-453-160 | Adson Bayonet | 1.0, 2.0 mm | 16.0 cm, 6 1/4" |
| TC20-456-191 | Adson Bayonet | 1.0, 2.0 mm | 19.1 cm, 7 1/2" |
| TC20-459-160 | Adson Bayonet (Ang. Up) | 1.0, 2.0 mm | 16.0 cm, 6 1/4" |
| TC20-462-191 | Adson Bayonet (Ang. Up) | 1.0, 2.0 mm | 19.1 cm, 7 1/2" |
| TC20-465-210 | Adson Bayonet (Ang. Up) | 1.0, 2.0 mm | 21.0 cm, 8 1/4" |
| TC20-468-241 | Remington Hobb | 1.0, 2.0 mm | 24.1 cm, 9 1/2" |

European Monopolar (Female 4.0 mm Connector)

| Cat. # | Description | Tip (mm) | Size |
|--------------|-------------------------|-------------|-----------------|
| TC20-320-120 | Adson | 1.0 mm | 12.0 cm, 4 1/4" |
| TC20-323-150 | McIndoe | 1.0, 2.0 mm | 15.0 cm, 6" |
| TC20-326-150 | Gillies (1 x 2) | -- | 15.0 cm, 6" |
| TC20-327-178 | Gillies (1 x 2) | -- | 17.8 cm, 7" |
| TC20-330-200 | Gillies (1 x 2) | -- | 20.0 cm, 8" |
| TC20-333-240 | Gillies (1 x 2) | -- | 24.0 cm, 9" |
| TC20-336-150 | Dissecting Forceps | 1.0, 2.0 mm | 15.0 cm, 6" |
| TC20-339-178 | Dissecting Forceps | 1.0, 2.0 mm | 17.8 cm, 7" |
| TC20-342-210 | Dissecting Forceps | 1.0, 2.0 mm | 21.0 cm, 8 1/4" |
| TC20-345-241 | Dissecting Forceps | 1.0, 2.0 mm | 24.1 cm, 9 1/2" |
| TC20-348-178 | Gerald (1 x 2) | -- | 17.8 cm, 7" |
| TC20-351-160 | Adson Bayonet | 1.0, 2.0 mm | 16.0 cm, 6 1/4" |
| TC20-354-191 | Adson Bayonet | 1.0, 2.0 mm | 19.1 cm, 7 1/2" |
| TC20-357-160 | Adson Bayonet (Ang. Up) | 1.0, 2.0 mm | 16.0 cm, 6 1/4" |
| TC20-360-191 | Adson Bayonet (Ang. Up) | 1.0, 2.0 mm | 19.1 cm, 7 1/2" |
| TC20-363-210 | Adson Bayonet (Ang. Up) | 1.0, 2.0 mm | 21.0 cm, 8 1/4" |
| TC20-366-241 | Remington Hobb | 1.0, 2.0 mm | 24.1 cm, 9 1/2" |

European Monopolar (Male 4.0 mm Connector)

| Cat. # | Description | Tip (mm) | Size |
|--------------|-------------------------|-------------|-----------------|
| TC20-369-120 | Adson | 1.0 mm | 12.0 cm, 4 1/4" |
| TC20-372-150 | McIndoe | 1.0, 2.0 mm | 15.0 cm, 6" |
| TC20-375-150 | Gillies (1 x 2) | -- | 15.0 cm, 6" |
| TC20-378-178 | Gillies (1 x 2) | -- | 17.8 cm, 7" |
| TC20-381-200 | Gillies (1 x 2) | -- | 20.0 cm, 8" |
| TC20-384-240 | Gillies (1 x 2) | -- | 24.0 cm, 9" |
| TC20-387-150 | Dissecting Forceps | 1.0, 2.0 mm | 15.0 cm, 6" |
| TC20-390-178 | Dissecting Forceps | 1.0, 2.0 mm | 17.8 cm, 7" |
| TC20-393-210 | Dissecting Forceps | 1.0, 2.0 mm | 21.0 cm, 8 1/4" |
| TC20-396-241 | Dissecting Forceps | 1.0, 2.0 mm | 24.1 cm, 9 1/2" |
| TC20-399-178 | Gerald (1 x 2) | -- | 17.8 cm, 7" |
| TC20-402-160 | Adson Bayonet | 1.0, 2.0 mm | 16.0 cm, 6 1/4" |
| TC20-405-191 | Adson Bayonet | 1.0, 2.0 mm | 19.1 cm, 7 1/2" |
| TC20-408-160 | Adson Bayonet (Ang. Up) | 1.0, 2.0 mm | 16.0 cm, 6 1/4" |
| TC20-411-191 | Adson Bayonet (Ang. Up) | 1.0, 2.0 mm | 19.1 cm, 7 1/2" |
| TC20-414-210 | Adson Bayonet (Ang. Up) | 1.0, 2.0 mm | 21.0 cm, 8 1/4" |
| TC20-417-241 | Remington Hobb | 1.0, 2.0 mm | 24.1 cm, 9 1/2" |



TUNER WARWICK
TC20-474-178
 17.8 cm, 7" (Cuv.)

SURGICAL SCISSORS
TC20-486-178
 17.8 cm, 7" (Cuv.)

WILSON HEY FORCEPS
TC20-498-178
 17.8 cm, 7" (Cuv.)

ARTERY FORCEPS
TC20-510-178
 17.8 cm, 7" (Cuv.)

TUNER WARWICK
TC20-471-178
 17.8 cm, 7" (Str.)

SURGICAL SCISSORS
TC20-483-178
 17.8 cm, 7" (Str.)

WILSON HEY FORCEPS
TC20-495-178
 17.8 cm, 7" (Str.)

ARTERY FORCEPS
TC 20-507-178
 17.8 cm, 7" (Str.)

RICHEE FORCEPS
TC20-519-178
 17.8 cm, 7"

English Monopolar
 (Male 4.8 mm
 Connector)


European monopolar
 (Female 4.0 mm
 Connector)

| Cat. # | Description | Size |
|--------------|--------------------------------|-------------|
| TC20-471-178 | Turner-Warwick Scissors (Str.) | 17.8 cm, 7" |
| TC20-474-178 | Turner-Warwick Scissors (Crv.) | 17.8 cm, 7" |
| TC20-477-230 | Turner-Warwick Scissors (Str.) | 23.0 cm, 9" |
| TC20-480-230 | Turner-Warwick Scissors (Crv.) | 23.0 cm, 9" |
| TC20-483-178 | Surgical Scissors (Str.) | 17.8 cm, 7" |
| TC20-486-178 | Surgical Scissors (Crv.) | 17.8 cm, 7" |
| TC20-489-230 | Surgical Scissors (Str.) | 23.0 cm, 9" |
| TC20-492-230 | Surgical Scissors (Crv.) | 23.0 cm, 9" |
| TC20-495-178 | Wilson Hey Forceps (Str.) | 17.8 cm, 7" |
| TC20-498-178 | Wilson Hey Forceps (Crv.) | 17.8 cm, 7" |
| TC20-501-230 | Wilson Hey Forceps (Str.) | 23.0 cm, 9" |
| TC20-504-230 | Wilson Hey Forceps (Crv.) | 23.0 cm, 9" |
| TC20-507-178 | Artery Forceps (Str.) | 17.8 cm, 7" |
| TC20-510-178 | Artery Forceps (Crv.) | 17.8 cm, 7" |
| TC20-513-230 | Artery Forceps (Str.) | 23.0 cm, 9" |
| TC20-516-230 | Artery Forceps (Crv.) | 23.0 cm, 9" |
| TC20-519-178 | Richie Forceps | 17.8 cm, 7" |
| TC20-522-200 | Richie Forceps | 20.0 cm, 8" |
| TC20-525-230 | Richie Forceps | 23.0 cm, 9" |

| Cat. # | Description | Size |
|---------------|--------------------------------|-------------|
| TC 20-528-178 | Turner-Warwick Scissors (Str.) | 17.8 cm, 7" |
| TC 20-531-178 | Turner-Warwick Scissors (Crv.) | 17.8 cm, 7" |
| TC 20-534-230 | Turner-Warwick Scissors (Str.) | 23.0 cm, 9" |
| TC 20-537-230 | Turner-Warwick Scissors (Crv.) | 23.0 cm, 9" |
| TC 20-540-178 | Surgical Scissors (Str.) | 17.8 cm, 7" |
| TC 20-543-178 | Surgical Scissors (Crv.) | 17.8 cm, 7" |
| TC 20-546-230 | Surgical Scissors (Str.) | 23.0 cm, 9" |
| TC 20-549-230 | Surgical Scissors (Crv.) | 23.0 cm, 9" |
| TC 20-552-178 | Wilson Hey Forceps (Str.) | 17.8 cm, 7" |
| TC 20-555-178 | Wilson Hey Forceps (Crv.) | 17.8 cm, 7" |
| TC 20-558-230 | Wilson Hey Forceps (Str.) | 23.0 cm, 9" |
| TC 20-561-230 | Wilson Hey Forceps (Crv.) | 23.0 cm, 9" |
| TC 20-564-178 | Artery Forceps (Str.) | 17.8 cm, 7" |
| TC 20-567-178 | Artery Forceps (Crv.) | 17.8 cm, 7" |
| TC 20-570-230 | Artery Forceps (Str.) | 23.0 cm, 9" |
| TC 20-573-230 | Artery Forceps (Crv.) | 23.0 cm, 9" |
| TC 20-576-178 | Richie Forceps | 17.8 cm, 7" |
| TC 20-579-200 | Richie Forceps | 20.0 cm, 8" |
| TC 20-582-230 | Richie Forceps | 23.0 cm, 9" |



HOLSTER
Reusable
TC20-606-000



TIP CLEANER
Reusable
TC20-609-000



FINGER-SWITCH PENCIL
Single Use
TC20-585-003 (3m cable)
TC20-588-005 (5m cable)



ROCKER-SWITCH PENCIL
Single Use
TC20-591-003 (3m cable)



FOOT CONTROL PENCIL
Single Use
TC20-594-003 (3m cable)
TC20-597-005 (5m cable)



ROCKER-SWITCH PENCIL
Reusable
TC20-600-003 (3m silicone cable)



ELECTRODE HANDLE
Reusable
TC20-603-003 (3m silicone cable)





**SCHROEDER
TANACULUM**
TC20-612-240
24 cm, 9 1/2"



POZZI
TC20-615-260
26 cm, 10 1/4"



**SINGLE TOOTH
TANACULUM (24 cm)**
TC20-618-240
24 cm, 9 1/2"



CAMPION FORCEPS
TC20-621-220
22 cm, 8 3/4"



RING FORCEPS
TC20-624-220
22 cm, 8 3/4"



KELLY
TC20-627-317
31.7 cm, 12 1/2"



GELPI RETRACTOR
TC20-630-001



RIGBY RETRACTOR
TC20-633-002



SCHROEDER TENACULUM
TC20-636-003



IRIS HOOK
TC20-639-004



EMMETT HOOK
TC20-642-005



STRAIGHT HOOK
TC20-645-006



TWO-PRONG HOOK
TC20-648-007



THREE-PRONG HOOK
TC20-651-008



**DRESSING
FORCEPS**
TC20-654-210
21 cm, 8 1/4"
TC20-657-250
25 cm, 9 1/4"

**TISSUE
FORCEPS**
TC20-660-210
21 cm, 8 1/4"
TC20-663-250
25 cm, 9 1/4"

IRIS HOOK
TC20-666-250
25 cm, 9 1/4"

EMMET HOOK
TC20-669-250
25 cm, 9 1/4"


STRAIGHT HOOK
TC20-672-250
25 cm, 9 1/4"

**TWO-PRONG
HOOK**
TC20-675-250
25 cm, 9 1/4"

**THREE-PRONG
HOOK**
TC20-678-250
25 cm, 9 1/4"



TISCHLER BIOPSY PUNCH FORCEPS
TC20-681-210
 21 cm, 8 1/4"



LATERAL VAGINAL RETRACTOR
TC20-684-210
 21 cm, 8 1/4"
 Leep / Lietz



LATERAL VAGINAL RETRACTOR
TC20-687-001



KOGAN ENDOCERVICAL SPECULUM
TC20-690-001



TC20-693-002
 3 mm



TC20-696-003
 4 mm



TC20-699-004
 5 mm



LATERAL RETRACTOR
 (wide open)
TC20-702-001



RIGBY RETRACTOR
TC20-705-002



GELPI RETRACTOR
TC20-708-003



COLLIN SPECULUM
With Smoke Evacuation Tube
TC20-711-001 Large
TC20-714-002 Medium
TC20-717-003 Small



CUSCO SPECULUM
With Smoke Evacuation Tube
TC20-720-001 Large
TC20-723-002 Medium
TC20-726-003 Small



GRAVES VAGINAL SPECULUM
TC20-729-001 SBC-Small
TC20-732-002 MBC-Medium
TC20-735-003 LBC-Large
TC20-738-004 BCX-Large



PEDERSON VAGINAL SPECULUM
TC20-741-001 SBC-Small
TC20-744-002 MBC-Medium
TC20-747-003 LBC-Large
TC20-750-004 BCX-Large



VAGINAL SPECULUM
TC20-765-001 Small
TC20-768-002 Medium
TC20-771-003 Large



WEISMAN SPECULUM
(Left Open)
TC20-753-001 Large
TC20-756-002 Medium



WEISMAN SPECULUM
(Right Open)
TC20-759-001 Large
TC20-762-002 Medium



VAGINAL SPECULUM
TC20-774-001 Small
TC20-777-002 Medium
TC20-780-003 Large



GRAVES VAGINAL SPECULUM
With Smoke Evacuation Tube
TC20-783-001 Large
TC20-786-002 Medium
TC20-789-003 Small



Large **TC20-783-001**
 Medium **TC20-786-002**
 Small **TC20-789-003**



PEDERSON VAGINAL SPECULUM
With Smoke Evacuation Tube
TC20-792-001 Large
TC20-795-002 Medium
TC20-798-003 Small



Large **TC20-792-001**
 Medium **TC20-795-002**
 Small **TC20-798-003**



VU-MAX SPECULUM
With Smoke Evacuation Tube
TC20-801-001 Graves (Medium)
TC20-804-002 Pederson (Medium)



TC20-801-001 Graves (M)
TC20-804-002 Pederson (M)



VU-MORE (Graves) SPECULUM
With Smoke Evacuation Tube
TC20-807-001 Large
TC20-810-002 Medium



TC20-807-001 Large
TC20-810-002 Medium



GRAVES VAGINAL SPECULUM
With Smoke Evacuation Tube
TC20-813-001 Large
TC20-816-002 Medium
TC20-819-003 Small



COLLIN SPECULUM
With Smoke Evacuation Tube
TC20-831-001 Large
TC20-834-002 Medium
TC20-837-003 Small



PEDERSON VAGINAL SPECULUM
With Smoke Evacuation Tube
TC20-822-001 Large
TC20-825-002 Medium
TC20-828-003 Small



LATERAL VAGINAL SPECULUM
TC20-840-001



KOGAN ENDOCERVICAL SPECULUM
TC20-843-001



TC20-846-002
3 mm








TC20-849-003
4 mm









TC20-852-004
5 mm








| | | | | | |
|------------|---|---|---|---|---|
| |  |  |  |  |  |
| Size | (1.0x0.4cm) | (1.0x0.5cm) | (1.0x0.7cm) | (1.0x0.8cm) | (1.0x1.0cm) |
| Disposable | TC20-855-001 | TC20-861-003 | TC20-867-005 | TC20-873-007 | TC20-879-009 |
| Reusable | TC20-858-002 | TC20-864-004 | TC20-870-006 | TC20-876-008 | TC20-882-010 |









| | | | | | | |
|------------|---|---|---|---|---|---|
| |  |  |  |  |  |  |
| Size | (1.5x0.5cm) | (1.5x0.7cm) | (1.5x0.8cm) | (1.5x1.0cm) | (1.5x1.2cm) | (1.5x1.5cm) |
| Disposable | TC20-885-001 | TC20-891-003 | TC20-897-005 | TC20-903-007 | TC20-909-009 | TC20-915-011 |
| Reusable | TC20-888-002 | TC20-894-004 | TC20-900-006 | TC20-906-008 | TC20-912-010 | TC20-918-012 |



| | | | | | |
|------------|---|---|---|---|---|
| |  |  |  |  |  |
| Size | (1.0x0.5cm) | (1.0x0.7cm) | (1.0x0.8cm) | (1.0x1.0cm) | (1.0x1.5cm) |
| Disposable | TC20-921-001 | TC20-927-003 | TC20-933-005 | TC20-939-007 | TC20-945-009 |
| Reusable | TC20-924-002 | TC20-930-004 | TC20-936-006 | TC20-942-008 | TC20-948-010 |



| | | | | | | |
|------------|---|---|---|---|---|---|
| |  |  |  |  |  |  |
| Size | (2.0x0.8cm) | (2.0x1.0cm) | (2.0x1.2cm) | (2.0x1.3cm) | (2.0x1.5cm) | (2.0x2.0cm) |
| Disposable | TC20-951-001 | TC20-957-003 | TC20-963-005 | TC20-969-007 | TC20-975-009 | TC20-981-011 |
| Reusable | TC20-954-002 | TC20-960-004 | TC20-966-006 | TC20-972-008 | TC20-978-010 | TC20-984-012 |



BLADE ELECTRODE
 Disposable: TC20-987-001
 Reusable: TC20-990-002



BLADE ELECTRODE (crv.)
 Disposable: TC20-993-003
 Reusable: TC20-996-004



NEEDLE ELECTRODE
 Disposable: TC20-999-005
 Reusable: TC20-1002-006



NEEDLE ELECTRODE (crv.)
 Disposable: TC20-1005-007
 Reusable: TC20-1008-008



BALL ELECTRODE (2mm)
 Disposable: TC20-1011-009
 Reusable: TC20-1014-010



BALL ELECTRODE CRV. (2mm)
 Disposable: TC20-1017-011
 Reusable: TC20-1020-012



BALL ELECTRODE (3mm)
 Disposable: TC20-1023-013
 Reusable: TC20-1026-014



BALL ELECTRODE CRV. (3mm)
 Disposable: TC20-1029-015
 Reusable: TC20-1032-016




BALL ELECTRODE (4mm)
 Disposable: TC20-1035-017
 Reusable: TC20-1038-018



BALL ELECTRODE CRV. (4mm)
 Disposable: TC20-1041-019
 Reusable: TC20-1044-020



BALL ELECTRODE (5mm)
 Disposable: TC20-1047-021
 Reusable: TC20-1050-022



BALL ELECTRODE CRV. (5mm)
 Disposable: TC20-1053-023
 Reusable: TC20-1056-024



BLADE ELECTRODE (Shaft 13.0cm)
 Disposable: TC20-1059-025
 Reusable: TC20-1062-026



NEEDLE ELECTRODE (Shaft 13.0cm)
 Disposable: TC20-1065-027
 Reusable: TC20-1068-028



BALL ELECTRODE 3mm (Shaft 13.0cm)
 Disposable: TC20-1071-029
 Reusable: TC20-1074-030



BALL ELECTRODE 5mm (Shaft 13.0cm)
 Disposable: TC20-1077-031
 Reusable: TC20-1080-032



CONIZATION ELECTRODE
 Disposable: TC20-1083-033
 Reusable: TC20-1086-034



CIRCLE LOOP ELECTRODE (4.0cm)
 Reusable: TC20-1089-035




CIRCLE LOOP ELECTRODE (6.0cm)
 Reusable: TC20-1092-036



CIRCLE LOOP ELECTRODE (8.0cm)
 Reusable: TC20-1095-037



DIAMOND ELECTRODE
 Reusable: TC20-1098-038



FINE NEEDLE ELECTRODE
 Reusable: TC20-1101-039



TC20-1104-001 - Silicone bipolar cable 10' (3.0m)
 Instrument Connector: USA 2 - Pin Plug
 Generator Connector: 4 mm Banana Plug
 Compatibility: Valleylab, Conmed, Codman, Olympus & Most Bipolar Generator




TC20-1125-008 - Silicone Bipolar Cable 10' (3.0m)
 Instrument Connector: USA 2 - Pin Plug
 Generator Connector: Pin Distance 29 mm
 Compatibility: Valleylab, EMC, Bolley



TC20-1113-004 - Silicone Bipolar Cable 10' (3.0m)
 Instrument Connector: USA 2 - Pin Plug
 Generator Connector: Martin Type
 Compatibility: Aesculap, Martin, Berchtold, Wolf



TC20-1110-003 - Silicone European Bipolar Cable 10' (3.0m)
 Instrument Connector: European Flat Plug
 Generator Connector: Martin Type
 Compatibility: Aesculap, Martin, Berchtold, Wolf



TC20-1122-007 - Silicone European Bipolar Cable 10' (3.0m)
 Instrument Connector: European Flat Plug
 Generator Connector: Pin Distance 29 mm
 Compatibility: Valleylab, EMC, Bolley



TC20-1119-006 - Silicone Bipolar Cable 10' (3.0m)
 Instrument Connector: USA 2 - Pin Plug
 Generator Connector: Erbe Type
 Compatibility: Erbe, Downs, Siemens




TC20-1107-002 - Silicone european bipolar cable 10' (3.0m)
 Instrument Connector: European Flat Plug
 Generator Connector: 4 mm Banana Plug
 Compatibility: Valleylab, Conmed, Codman, Olympus & Most Bipolar Generator



TC20-1128-009 - Silicone Bipolar Cable 10' (3.0m)
 Instrument Connector: USA 2 - Pin Plug
 Generator Connector: Pin Distance 29 mm
 Compatibility: Valleylab, EMC, Bolley



TC20-1116-005 - Silicone Bipolar Cable 10' (3.0m)
 Instrument Connector: European Flat Plug
 Generator Connector: Erbe Type
 Compatibility: Erbe, Downs, Siemens



TC20-1131-001 - Silicone Monopolar cable 10' (3.0m)
 Instrument Connector: Screw Type
 Generator Connector: Standar Bovie 8 mm Plug
 Compatibility: Valleylab, Conmed, Bovie



TC20-1143-005 - Silicone Diathermy cable 10' (3.0m)
 Instrument Connector: 4.8 mm Female Fitting (English Diathermy)
 Generator Connector: Standar Bovie 8 mm Banana Plug
 Compatibility: Valleylab, Conmed, Bovie



TC20-1134-002 - Silicone Monopolar Endoscopic cable 10' (3.0m)
 Instrument Connector: 4 mm Female Fitting
 Generator Connector: Standar Bovie 8 mm Plug
 Compatibility: Valleylab, Conmed, Bovie



TC20-1140-004 - Silicone Diathermy cable 10' (3.0m)
 Instrument Connector: 4 mm Male Fitting (European Diathermy)
 Generator Connector: Standar 4 mm Banana Plug
 Compatibility: Most Es Generators



TC20-1137-003 - Silicone Monopolar Endoscopic cable 10' (3.0m)
 Instrument Connector: 4 mm Female Fitting
 Generator Connector: Standar 4 mm Banana Plug
 Compatibility: Most Es Generators




TC20-1146-006 - Silicone Diathermy cable 10' (3.0m)
 Instrument Connector: 4.8 mm Female Fitting (English Diathermy)
 Generator Connector: Standar 4 mm Banana Plug
 Compatibility: Most Es Generators

TTC
MEDICLINE



TTC
MEDICLINE



| | | | | | | | | |
|---|---|---|---|---|---|---|---|--|
| HOCKEY TC20-1149-120 12 cm, 4 3/4" | HOCKEY TC20-1152-120 12 cm, 4 3/4" | HOOK TC20-1155-180 18 cm, 7 1/4" | HOOK TC20-1158-180 18 cm, 7 1/4" | DAGGER TC20-1161-180 18 cm, 7 1/4" | DAGGER TC20-1164-120 12 cm, 4 3/4" | FRAZIER SUCTION TUBE 11 FRENCH TC20-1167-135 13.5 cm, 5 1/4" | FRAZIER SUCTION TUBE 11 FRENCH TC20-1170-205 20.5 cm, 8" | FRAZIER SUCTION TUBE 11 FRENCH TC20-1173-406 40.6 cm, 16" |
|---|---|---|---|---|---|---|---|--|

WITH SUCTION



L-HOOK ELECTRODE
TC20-1176-360
36 cm, 14 "



J-HOOK ELECTRODE
TC20-1179-360
36 cm, 14 "



NEEDLE ELECTRODE
TC20-1182-360
36 cm, 14 "



SPATULA ELECTRODE
TC20-1185-360
36 cm, 14 "



BALL ELECTRODE
TC20-1188-360
36 cm, 14 "

WITHOUT SUCTION



L-HOOK ELECTRODE
TC20-1191-360
36 cm, 14 "



J-HOOK ELECTRODE
TC20-1194-360
36 cm, 14 "



NEEDLE ELECTRODE
TC20-1197-360
36 cm, 14 "



SPATULA ELECTRODE
TC20-1200-360
36 cm, 14 "



BALL ELECTRODE
TC20-1203-360
36 cm, 14 "



Wazirbad Road, Ugoki, Sialkot-51310 PAKISTAN

Phone: +92-52-3561347-3553588

Phone: +92-52-3561202

Cell: +92-300-6139615

Cell: +92-322-6139615

Email: info@tcmedicline.com

Email: exporttcmedicline@gmail.com